2526
GALLONI

The future belongs to those who have a great story

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PRESSOVAC SERIES ACCESSORIES

- External high VACUUM PUMP.
- WATER PUMP supplied complete with water recirculating cooling tank.

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The system used by our Pressovac series casting machines was specially developed for dental labs where casting perfection is absolute: casting pieces have a very smooth surface, a very dense molecular structure and without defects like improper filling, porosity and shrinkage.

Casting cycle is fully automatic. The casting chamber is fully evacuated by vacuum before melting, followed by helium gas supply to create an inert atmosphere.

Helium gas has an excellent gas fluidity which provides a cleaning effect of the molten metal and mould and guarantees a smooth filling of metal even of the most intricate patterns.

The casting temperature is attained and low frequency pulses are given to vibrate and mix the molten metal by keeping it homogeneous. Temperature is controlled by an optical pyrometer up to 2100°C (3812°F) for platinum and steel and by an immersion type "K" thermocouple inside the graphite crucible up to 1150°C (2102°F) for gold, silver and bronze.

At this stage, just by pushing a button, the whole crucible/flask system will automatically rotate 90° and the metal, once poured in a homogeneous and controlled way into the flask, it is automatically pressurized to 3 bars by argon gas. This high and quick pressurization is indispensable to obtain a proper and homogeneous filling, even at low metal/flask temperatures, for gold, silver and bronze alloys and, above all, for platinum and special steels which have a very narrow melting interval, thus a very quick solidification. This results in a very smooth surface due to the effect of the vacuum and a dense, well compacted casting due to the "over pressure" of 3 bars, which benefits both thick and thin sections.

Features of the pressovac series

- Extremely robust water-cooled casting chamber
- Automatic tilt
- Rapid pressurization up to 3 bar, allowing lower metal and flask temperatures
- Automatic lock of cover

Competitors

Most of our competitors' machines have weak casting chambers which cannot withstand 3 bar pressure supplied in a very short time, otherwise the cover will blow out. Therefore gas pressure has to be applied slowly to the molten metal after pouring. This requires a significant temperature increment of both flask and metal to achieve adequate filling, resulting in porosity and rough surface of the casted piece. A further disadvantage of competitors' machines is the manual tilting, an action affected by the operator skill: incomplete filling if the pouring is slow, metal spillage if too rapid.

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The "DUAL PYROMETER" version of our well-proven Pressovac machine features:

- optical infrared pyrometer for temperatures up to 2100°C (3812°F) for steel and platinum.
- immersion "K" thermocouple for temperatures up to 1150°C (2102°F) for gold, silver, brass, bronze.

Example of ministructure

The ministructure consists of the casting machine, the burnout furnace, the investment mixer and the water jet cabin; it requires a small area on a bench top. Goldsmiths, jewellery designers and studios can cast by themselves high quality intricate patterns in a short time. An ideal way to produce from rapid-prototyping.

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
Touch-sceen capable to store 50 casting processes and to display the actual parameters of pressure, temperature and power.

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EXAMPLE OF FILIGREE CASTING FROM RAPID PROTOTYPING MODELS (courtesy by DWS - Italy)

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Equipped with an "universal cradle" for any flask diameter (between 70 mm ø minimum and 90 mm ø maximum) by 90 mm height.

PRESSOVAC DUAL PYROMETER TECHNICAL FEATURES / CARATTERISTICHE TECNICHE

ELECTRICAL / ALIMENTAZIONE	singlephase 230 V 50/60 Hz - 3.5 kW
FLASKS / CILINDRI	from 70 to 90 mm ø, up to 90 mm h
COOLING LIQUID / RAFFREDDAMENTO	2 l/minute 3~4 bar with closed circuit water pump
COMPRESSED AIR / ARIA COMPRESSA	6~7 bar
GAS / GAS	Helium, Argon, Nitrogen
EXTERNAL VACUUM PUMP / POMPA VUOTO ESTERNA	8 m ³ /h - 0.5 mbar
OVERAIL DIMENSIONS AND WEIGHT / DIMENSIONE PESO	600 x 500 mm x 580 mm h.- 90 kg

In compliance with CE regulation - Warranty 12 months only if Galloni original consumables are used

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
The "MAX" version of our Pressovac machine is free-standing and offers the same proven features and functions as the standard Pressovac but accepting a larger flask and with an enhanced casting capacity (flask size up to 100 mm ø x 150 mm h.)

All functions are touch-screen controlled.

The machine features:

- Optical infrared pyrometer for temperatures up to 2100°C (3812°F) for steel and platinum
- Immersion "K" thermocouple for temperatures up to 1150°C (2102°F) for gold, silver, brass, bronze.
- Option: immersion "S" thermocouple for temperatures up to 1450°C (2642°F) for gold alloys with palladium.

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ELECTRICAL / ALIMENTAZIONE	threephase 230/400/415 V 50/60 Hz – 6 kW
FLASKS / CILINDRI	from 70 to 100 mm ø, up to 150 mm h
COOLING LIQUID / RAFFREDDAMENTO	3 l/minute 3-4 bar with closed circuit water pump
COMPRESSED AIR / ARIA COMPRESSA	6~7 bar
GAS / GAS	Helium, Argon, Nitrogen
EXTERNAL VACUUM PUMP / POMPA VUOTO ESTERNA	8 m ³ /h – 0.5 mbar
OVERALL DIMENSIONS AND WEIGHT / DIMENSIONE PESO	680 x 490 mm x 1200 mm h. – 160 kg

In compliance with CE regulation – Warranty 12 months only if Galloni original consumables are used

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PRESSOVAC MAX XL is the high capacity version of the series PRESSOVAC.

It accepts flask sizes up to 150 mm \varnothing and 250 mm h and is controlled by a colour touch screen.

The machine features:

- Optical infrared pyrometer for temperatures up to 2100°C (3812°F) for steel and platinum
- Immersion "S" thermocouple for temperatures up to 1150°C (2102°F) for gold, silver, brass, bronze.
- Option: immersion "S" thermocouple for temperatures up to 1450°C (2642°F) for gold alloys with palladium.

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ELECTRICAL / ALIMENTAZIONE	threephase 230/400/415 V-50/60 Hz -10 kW
FLASKS / CILINDRI	from 100 to 150 mm \varnothing , from 120 to 250 mm h
COOLING LIQUID / RAFFREDDAMENTO	water 3 l/minute at 3 bar with closed circuit water pump
COMPRESSED AIR / ARIA COMPRESSA	7 bar
GAS / GAS	Helium, Argon, Nitrogen
EXTERNAL VACUUM PUMP / POMPA VUOTO ESTERNA	60 m ³ /h - 0.5 mbar
OVERALL DIMENSIONS AND WEIGHT / DIMENSIONE PESO	1000 x 700 mm x h 1250 mm - 330 kg

In compliance with CE regulation - Warranty 12 months only if Galloni original consumables are used

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PRESSOVAC DUAL PYROMETER

Au KIT

A100.566	Graphite crucible cemented with its container for gold and silver (20 cc)
A100.509	Lid cover for crucible
A100.580	Graphite thermocouple sheath
A100.056	Spare thermocouple type "K" for 1150°C (2102°F)
P06101	Melting coil (101)
P06099	Crucible holding plate complete with crucible lid cover holder mechanism

Pt/Pd/Steel KIT

A100.559	Ceramic crucible for platinum (10 cc)
A100.564	Silicon carbide crucible (20 cc)
P06102	Melting coil (102)
P06098	Crucible holding plate complete with crucible lid cover holder mechanism

PRESSOVAC MAX

Au KIT

A100.594	Graphite crucible cemented with its container for gold and silver (60 cc)
A100.577	Lid cover for crucible
A100.575	Graphite thermocouple sheath
A083.055	Spare thermocouple type "K" for 1150°C (2102°F)
P06110	Melting coil (103)
P06113	Crucible holding plate complete with crucible lid cover holder mechanism

Pt KIT

A100.512	Ceramic crucible for platinum (20 cc)
P06111	Melting coil (104)
P06114	Crucible holding plate complete with crucible lid cover holder mechanism

Steel/Pd KIT

A100.568	Silicon carbide crucible (40 cc)
P06111	Melting coil (104)
P06115	Crucible holding plate complete with crucible lid cover holder mechanism

Steel/Pd KIT (Big)

A100.578	Silicon carbide crucible (75 cc)
P06112	Melting coil (105)
P06116	Crucible holding plate complete with crucible lid cover holder mechanism

PRESSOVAC MAX XL

Au KIT

A100.434	Graphite crucible cemented with its container for gold and silver (60 cc)
A100.410	Lid cover for crucible
	Graphite thermocouple sheath
	Spare thermocouple type "K" for 1150°C (2102°F)
P06071	Melting coil (3668)
P06122	Crucible holding plate complete with crucible lid cover holder mechanism

Pt KIT

A100.649	Ceramic crucible for platinum
P06117	Melting coil (106)
P06121	Crucible holding plate complete with crucible lid cover holder mechanism

Pd/Steel KIT

A100.626	Silicon carbide crucible (60 cc)
P06071	Melting coil (3668)
P06120	Crucible holding plate complete with crucible lid cover holder mechanism

Since 1957
The technology in the
investment casting

ASEG
GALLONI

...there's no better way to cast!

ASEG GALLONI S.P.A.
Via Don Luigi Sturzo, 20 - 20078
San Colombano al Lambro (MI) - Italy
T. +39 0371 200233 / F. +39 0371 898705
info@galloni-aseg.com
www.galloni-aseg.com